

Product Information

Polyacetal agent, Chondroitin sulfate-polyacetal, Purity $\geq 95\%$

Cat. No.: X25-05-YM023

Size: 100 mg; 250 mg; 500 mg; 1 g; 5 g

Synonym: Polyacetal-chondroitin sulfate; Polyacetal-chondroitin sulfate; Chondroitin sulfate-polyacetal

This product is for research use only and is not intended for diagnostic use.

Product Information

Description	Polyacetal-chondroitin sulfate graft copolymer is engineered through acetalization reactions, creating acid-degradable bonds between synthetic polymer chains and chondroitin sulfate moieties.
Glycan Name	Chondroitin sulfate
Glycan Structure	The glycan structure of chondroitin sulfate (CS) is a sulfated glycosaminoglycan (GAG) composed of repeating disaccharide units. Each unit consists of: N-acetyl-D-galactosamine (GalNAc) ($\beta 1 \rightarrow 4$ linked) D-glucuronic acid (GlcA) ($\beta 1 \rightarrow 25$ linked)
Source	Chemical synthesis
Functional Group	Polyacetal
Form	Solid or powder
Purity	$\geq 95\%$
Impurities	No visible impurities to the naked eye.
Solubility	This product is soluble in most organic solvents, such as DCM, DMF, DMSO, and THF, and exhibits excellent solubility in water.
Identity	Confirmed by NMR.
Stability	It is stable under its storage temperature.
Quality Level	Research level
Applications	Chondroitin sulfate-polyacetal can be used for its potential to develop acid-degradable polymers for lysosomal drug release.
Storage	Store at -20°C , protect from light and moisture.