

Product Information

FITC agent, Heparin-FITC, Purity $\geq 95\%$

Cat. No.: X25-11-YM103

Size: 100 mg; 250 mg; 500 mg; 1 g; 5 g

Synonym: FITC-PEG-heparin

This product is for research use only and is not intended for diagnostic use.

Product Information

Description	A fluorescently labeled heparin conjugate has fluorescein isothiocyanate (FITC) covalently attached, facilitating visualization and tracking of heparin distribution in cellular and tissue imaging studies.
Functional Group	FITC
Form	Solid or powder
Purity	$\geq 95\%$
Impurities	No visible impurities to the naked eye.
Solubility	This product is soluble in most organic solvents, such as DCM, DMF, DMSO, and THF, and exhibits excellent solubility in water.
Identity	Confirmed by NMR.
Stability	It is stable under its storage temperature.
Quality Level	Research level
Applications	Heparin-FITC can enable fluorescent visualization of heparin localization and uptake in cellular studies.
Storage	Store at -20°C , protect from light and moisture.