



Product Information


4-O-Sulfate-GlcNAc β -Sp2-PAA-FITC

Cat. No.: X22-12-ZQ751

Size: 400 μ g; 800 μ g; 1 mg; 4 mg

Synonym: D-GlcNAc-4-O-sulfate-PAA-FITC; 4-suGlcNAc-C2-PAA-fluo

This product is for research use only and is not intended for diagnostic use.



Product Information

Molecular Weight	20 kDa
Glycan Structure	4-O-Sulfate-GlcNAc
Glycan MW	344.4 Da
Type	Polymeric
Label	FITC
Spacer	Sp2 (OCH ₂ CH ₂ NH)
Carrier	PAA
Form	Lyophilized powder
Purity	>95%, determined by HPLC and ¹ H NMR.
Applications	4-O-Sulfate-GlcNAc β -Sp2-PAA-FITC is a useful polymer for studying carbohydrate-protein interactions. This product can be used in flow cytometry (FC) and fluorescence microscopy.
Storage	Store at -20°C.

Safety Information

Personal Protective Equipment	Gloves, eyeshields, N95 mask
WGK Germany	WGK 3