

## Product Information


### Gal $\beta$ (1-4)[Gal $\beta$ (1-3)]GlcNAc $\beta$ -Sp3-PAA-FITC

**Cat. No.:** X22-12-ZQ815

**Size:** 400  $\mu$ g; 800  $\mu$ g; 1 mg; 4 mg

**Synonym:** (Gal)<sub>2</sub>-3,4GN-PAA-FITC; (Gal)<sub>2</sub>-3,4GN-C3-PAA-fluo

**This product is for research use only and is not intended for diagnostic use.**



#### Product Information

<b>Molecular Weight</b>	20 kDa
<b>Glycan Structure</b>	Gal $\beta$ (1-4)[Gal $\beta$ (1-3)]GlcNAc
<b>Glycan MW</b>	602.6 Da
<b>Type</b>	Polymeric
<b>Label</b>	FITC
<b>Spacer</b>	Sp3 (OCH <sub>2</sub> CH <sub>2</sub> CH <sub>2</sub> NH)
<b>Carrier</b>	PAA
<b>Form</b>	Lyophilized powder
<b>Purity</b>	>95%, determined by HPLC and <sup>1</sup> H NMR.
<b>Applications</b>	Gal $\beta$ (1-4)[Gal $\beta$ (1-3)]GlcNAc $\beta$ -Sp3-PAA-FITC is a useful polymer for studying carbohydrate-protein interactions. This product can be used in flow cytometry (FC) and fluorescence microscopy.
<b>Storage</b>	Store at -20°C.

#### Safety Information

<b>Personal Protective Equipment</b>	Gloves, eyeshields, N95 mask
<b>WGK Germany</b>	WGK 3