



Product Information

Schiff's fuchsin-sulfite reagent

Cat. No.: XGB977

Size: 500 mL; 1 L

NACRES: NA.47

This product is for research use only and is not intended for diagnostic use.

Product Information

Molecular Weight	383.46
Empirical Formula	C ₂₀ H ₂₁ N ₃ O ₃ S
Form	liquid
Purity	suitable for detection of glycoproteins
Quality Level	100
Applications	Schiff's fuchsin-sulfite reagent can be used: • to stain ovaries and testis • to stain undifferentiated sarcomas (USARCs) • to stain tomato and Medicago primary roots
Storage	2-8°C

Safety Information

Personal Protective Equipment	Eyeshields, Gloves, half-mask respirator (US), multi-purpose combination respirator cartridge (US)
WGK Germany	WGK 3
Hazard Statements	H290 - H350